

The Seprion ligand only binds aggregated proteins. The Seprion ligand is a unique tool for the detection, investigation and study of protein aggregates and Protein Aggregation Diseases.

Seprion application areas

Investigation and study of disease

Alzheimer's Disease
 Parkinson's Disease
 Huntington's Disease
 Prion disease (CJD, BSE, scrapie etc.)
 And many more.....

Aggregated proteins studied

β -amyloid
 tau
 α -synuclein
 Huntingtin
 Prion protein
 List continues to grow

Application to studies *in vitro*

Cell culture studies of disease eg drug development

Application to studies *in vivo*

Animal models of disease eg drug efficacy and clinical studies

Investigation of the design, manufacture and stability of therapeutic proteins

Therapeutic proteins eg monoclonal antibodies

Investigation of aggregated proteins in cancer

For example, study of aggregated p53

Seprion is an established technology


Established patent protected technology with many publications in leading scientific journals eg

J.P. Clewley et al. Prevalence of disease related prion protein in anonymous tonsil specimens in Britain: cross sectional opportunistic survey. 2009. *BMJ* 338:b1442

K. Sathasivam et al. Identical oligomeric and fibrillar structures captured from the brains of R6/2 and knock-in mouse models of Huntington's disease. (2010); *Human Mol. Gen.* 19(1):65-78

Used by University of London, Lancaster University, Veterinary Laboratory Agency and many more researchers for research and drug screening. Used commercially by IDEXX Laboratories Inc


Seprion can be provided coated onto magnetic beads or onto microplates




Selected examples of research using Seprion

A. Samples from two normal and three confirmed Alzheimer’s Disease brains were homogenised and incubated with Seprion-coated magnetic beads.

After washing, captured proteins were eluted and analysed for abnormal 1-40 and 1-42 β -amyloid and tau by commercial ELISA




and 1-42 β -amyloid and abnormal tau protein compared to the brains from age matched controls without Alzheimer’s Disease



The Seprion-coated beads demonstrated that the brains of the Alzheimer’s Disease patients contained large amounts of both abnormal 1-40

B. Brains of Alzheimer’s Disease model mice at various ages were homogenised. The abnormal β -amyloid protein was separated from the normal protein using the Seprion-coated beads and investigated by ELISA.

C. Seprion-coated beads (PAD beads) were used to capture brain homogenates from wild type (WT) and transgenic (R6/2) mice expressing the human huntingtin gene. Captured proteins were analysed by Western blotting with an anti-huntingtin antibody.



Homogenate from a transgenic mouse brain that was not captured by PAD beads was run as a control.

D. Kits based on Seprion received USDA approval for BSE and CWD testing, EU approval for BSE and scrapie testing and have been used to screen many millions of cows for BSE.



E. CJD Kits for human testing have been used by the Health Protection Agency for tonsil screening for vCJD.